en	ENGLISH	pg. 1	fr	FRANCAIS	pg. 12	sk	SLOVENČINA	pg. 24
ja	日本語	pg. 3	it	ITALIANO	pg. 14	sv	SVENSKA	pg. 26
cs	ČEŠTINA	pg. 5	ko	한국어	pg. 16			
da	DANSK	pg. 7	lv	LATVIEŠU	pg. 18			
de	DEUTSCH	pg. 8	nl	NEDERLANDS	pg. 20			
es	ESPAÑOL	pg. 10	no	NORSK	pg. 22			

WARNING

After use, dispose of the product and packaging in accordance with hospital and/or local government policy.

WARNING

For single use only. Do not reuse, reprocess, or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity and result in poor image quality or patient injury, illness, or death. Reuse, reprocessing, or resterilization may also contaminate the device and result in patient infection which may lead to patient illness or death. Infraredx makes no performance claim for product that is reused, reprocessed, or resterilized.

Materials and Equipment (*Sold Separately)

Dualpro IVUS+NIRS Imaging Catheter
Makoto Controller Sterile Barrier
Priming Accessories
*Makoto IVUS+NIRS Imaging System
*Makoto IVUS+NIRS Imaging System User Manual
*Pre-formed guide catheter (≥ 0.067 inches ≥ 1.70 mm I.D. min., 6F)
*Guidewire, 0.014 inches (0.36mm) maximum diameter
*Heparinized normal saline

Instructions for Use

Refer to the Makoto IVUS+NIRS Imaging System User Manual for guidance on preparation of the Makoto IVUS+NIRS Imaging System and Makoto Controller.

NOTE:

Information related to the safety and clinical performance (SSCP) of the Makoto™ Intravascular Imaging System is available in the European database on medical devices (Eudamed), where it is linked to the Basic UDI-DI located beneath the bar code next to (01) on the product label.

Eudamed Website: <https://ec.europa.eu/tools/eudamed>

Inspection Prior to Use

Prior to use, ensure that the catheter and sterile package have not been damaged or breached and that no particulate is present inside the catheter. Prior to scanning, all equipment to be used during the procedure should be carefully examined for defect to ensure proper performance.

WARNING

If Catheter Sterile Packaging is damaged or unintentionally opened before use, immediately dispose of the Packaging and its Contents and replace with a new Catheter Sterile Package.

Catheter Preparation

1. Using sterile technique, open the pouch and remove catheter from packaging. DO NOT remove catheter from the hoop until the hydrophilic coating has been activated;
2. Remove priming syringes from packaging;
3. Fill 10mL and 3mL syringes with heparinized saline and connect to the three-way stopcock.
4. Using the gray catheter cap, gently withdraw the catheter imaging core from the catheter handle until it stops.
5. Prime the Dualpro catheter TWICE using the 3 mL syringe. Use the 10 mL syringe as a reservoir to refill the 3 mL syringe. Priming the catheter in the hoop will activate the hydrophilic coating.
6. Using the gray cap, turn the colored side of the rotating hub to match the corresponding colored fixed handle button. Push the rotating and fixed handle together. An audible click will signal complete closure.

The catheter is now primed. The gray protective cap should remain on the catheter until ready to connect to the Makoto Controller.

CAUTION

DO NOT use any type of contrast media either in replacement of or in combination with the saline as priming medium.

Connecting the Dualpro™ Imaging Catheter to the Makoto™ Controller

WARNING

If Controller Sterile Barrier is contaminated or damaged at any time, immediately replace with a new Controller Sterile Barrier.

The catheter must be connected to the Makoto Controller by a sterile operator in the sterile field. However, before connecting the catheter, you must:

1. Prepare the Controller for sterile use, by covering the Makoto Controller with the sterile barrier (Refer to Makoto IVUS+NIRS Imaging System User Manual).
2. Confirm that the sterile barrier gasket is engaged to the Makoto Controller connection socket.
3. Remove the blue covering from the sterile barrier that covers the catheter connection socket by separating along the perforated tape.

CAUTION

Do not contaminate the fiber faces of the catheter or Makoto Controller during the connection process.

4. Remove the Dualpro Imaging Catheter from the hoop.
5. Remove the protective cap from the catheter by holding the fixed hub and pulling the gray cap.
6. Align the catheter, oriented with the colored side up, with the Makoto Controller catheter connection socket.
7. Insert the catheter into the socket and turn the catheter handle clockwise until the tactile response is felt.
8. Confirm the connection by observing the READY state indication on the Makoto Controller.
 WARNING
When connected to the Makoto IVUS+NIRS Imaging System, laser radiation is emitted from the distal end of the catheter. Do not stare into beam or view directly with optical instruments
9. With the connected catheter in the READY position, prime the system with heparinized saline using the 3 mL priming syringe.
10. Assess the IVUS image on the screen by enabling Live IVUS. Concentric bright rings will appear in the transverse IVUS image with a successful prime.

NOTE:

If the IVUS image is dark, prime with additional saline.

Clinical Benefit

The clinical benefit provided by the Makoto™ Intravascular Imaging System is in the identification of the lipid core plaques (LCPs). By identifying LCP, the Makoto™ Intravascular Imaging System provides information that is not as easily available from other technologies. The potential benefits of NIRS-IVUS-guided management are: 1) determination of the length of artery to stent, 2) selection of intensity of medical management, 3) identification and mitigation of the risk of embolic peri-stenting infarction, and 4) identification of patients and plaques at increased risk of future major adverse cardiac events (MACE).

Complications

The following complications may occur as a consequence of intravascular examination:

- Arterial dissection, injury, or perforation
- Acute myocardial infarction
- Ventricular fibrillation
- Total occlusion
- Unstable angina
- Air embolism
- Abrupt closure
- Death
- Cardiac Tamponade
- Thrombus Formation

Safety

Serious incidents involving the Makoto™ Intravascular Imaging System should immediately be reported to Infraredx and the Competent Authority of the Member State in which the user and/or patient is established. Infraredx or your local service provider contact information may be found on the last page.

- When satisfied with the catheter prime and image, press the STOP button on the Makoto Controller to halt motion.
- Prime the system with heparinized saline using the 3 mL priming syringe.

You are now ready to perform a scan with the Dualpro IVUS+NIRS Catheter.

CAUTION
Retracting the imaging core of a primed catheter outside the body may introduce air into the catheter sheath. Flush the catheter using the 3 mL priming syringe then advance the catheter imaging core to the fully distal position using the distal motion controls.

Introduce Catheter into Guide Catheter and Artery

WARNING
Never advance the Dualpro catheter without guidewire support.
Never advance the Dualpro catheter sheath without the imaging core advanced to its most distal or READY position.
Never advance or withdraw the Dualpro catheter without direct fluoroscopic visualization.
If resistance is encountered anytime during Dualpro imaging catheter sheath positioning, DO NOT pull, push, or rotate with excessive force.
Never advance the distal tip of the Dualpro catheter near the unsupportive end of the guidewire due to the risk of guidewire entanglement.

CAUTION
Guidewires that supply more stiffness near the distal tip are recommended.
An excessively tightened hemostasis valve may distort the IVUS image or cause permanent damage to the imaging core during rotation.

NOTE:

Tracking performance of the catheter is optimized when the imaging core is advanced to the READY position, the most distal position.

- Confirm that the catheter imaging core is in the READY position.
- Load the guidewire into the guidewire provision at the distal end of the catheter.
- Advance the imaging catheter into the guide catheter, up to the appropriate insertion depth marker on the catheter shaft. Tighten the hemostasis valve on the guide catheter. Tighten only enough to prevent fluid / blood leakage.
- Using fluoroscopy guidance, advance the Dualpro catheter over the guidewire until the Image Start position on the imaging core is beyond the region of interest.
- Use the "Live IVUS" or "Pullback" imaging controls on the Makoto Controller to perform either manual or automated imaging.

NOTE:

When starting from the READY position, the imaging tip will move within the sheath proximal to the TRS, approximately 5mm, before performing translated Live IVUS imaging or automated pullback imaging.

NOTE:

When Pullback or translated Live IVUS begins proximal to the READY position, there is no translation prior to imaging. Imaging will begin at the current location of the imaging tip.

If additional scans are intended:

- Advance the catheter imaging core to the READY position using the Image Core Positioning controls or the Return-to-READY control.
- Reposition the Dualpro Catheter on the guidewire as described in step 4.

Catheter Imaging

The catheter sheath, guide catheter, and guidewire must remain fixed when imaging. Refer to the Makoto IVUS+NIRS Imaging System User Manual for instructions on proper use of the Makoto IVUS+NIRS Imaging System Console and Controller for imaging.

When imaging is complete, advance the imaging core to the READY position, maintain the position of the guidewire and withdraw the catheter. If the imaging catheter is to be re-inserted, then prime and confirm the imaging core is in the READY position. To protect the catheter from damage between uses, coil the catheter and set aside until reuse. When ready to re-insert the catheter, prepare the catheter as previously instructed.

CAUTION
Inspect the guidewire exit port and imaging window prior to re-insertion to verify that no damage has occurred during withdrawal.

Catheter Disconnection

- Use the imaging core positioning or the Return-to-READY controls to advance the imaging core to the READY position.
- Turn the catheter fixed hub counter-clockwise until the colored button of the catheter is facing up.
- Remove the catheter fixed hub from the Makoto Controller.

Additional WARNINGS and CAUTIONS

WARNING
Do not advance the imaging core into a kinked sheath.

Care should be taken when a guidewire is exposed in a stented vessel. Catheters that do not encapsulate the guidewire may engage the stent between the junction of the catheter and guidewire.

Care should be taken when re-advancing a guidewire after stent deployment. A guidewire may exit between stent struts when re-crossing a stent that is not fully apposed. Subsequent advancement of the Dualpro catheter could cause entanglement between the catheter and the mal-apposed stent.

If resistance is met upon withdrawal of the catheter, verify source of resistance using fluoroscopy and ensure the catheter is not entangled in a stent or other interventional device then apply an appropriate catheter removal strategy.

If a Dualpro catheter sheath breach occurs during the procedure, do not advance the catheter core. Immediately remove the entire system of Dualpro, guide catheter, and guidewire using fluoroscopic guidance.

CAUTION
Do not kink or sharply bend the catheter at any time. This can cause drive cable failure. An insertion angle greater than 45 degrees is considered excessive.

Warranty

Seller hereby warrants that all goods sold by Seller are free from defects in material and workmanship. Seller warrants that under normal circumstances all disposable goods will remain free from defects for the shorter duration of their specified shelf life, or as required legally by regulatory authority with jurisdiction in the Buyer's place of business. Seller will replace, at its discretion, all defective disposable goods during this warranty period. After the termination of the warranty period, replacement goods may be purchased from Seller at its published rates. This sub-section states Seller's entire liability for Product warranty claims.

The above warranties shall not apply to defects or damage resulting from misuse, negligence or accident, including, without limitation: operation by untrained operators; operation outside of the environmental specifications for the disposable goods; and making of modifications to disposable goods not authorized by Seller. This warranty is not transferable without the express written consent of Seller. Seller is under no obligation to provide any services not specified in this Agreement EXCEPT AS EXPRESSLY SET FORTH IN THIS AGREEMENT, NO WARRANTY, WHETHER EXPRESSED OR IMPLIED, IS MADE WITH RESPECT TO THE DISPOSABLES INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ALL OF WHICH ARE EXPRESSLY DISCLAIMED BY SELLER.

Limitation of Liability

SELLER SHALL NOT HAVE LIABILITY BEYOND THE REMEDIES SET FORTH HEREIN, INCLUDING BUT NOT LIMITED TO ANY LIABILITY FOR THE EQUIPMENT NOT BEING AVAILABLE FOR USE, MISDIAGNOSIS, LOST PROFITS, LOSS OF BUSINESS OR FOR LOST OR CORRUPTED DATA. SELLER WILL NOT BE LIABLE FOR ANY CONSEQUENTIAL, SPECIAL, INDIRECT, OR PUNITIVE DAMAGES, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGES, OR FOR ANY CLAIM BY ANY THIRD PARTY. IN NO EVENT SHALL SELLER BE LIABLE OR RESPONSIBLE FOR ANY AMOUNT OF DAMAGES ABOVE THE AMOUNT INVOICED FOR THE GOODS WITH RESPECT TO WHICH DAMAGES ARE SOUGHT.

Claims-notice of Defects

In no event shall Seller be liable on any claims unless written notice of the claim is received by Seller within the following limitations: for shortages in weight and count, 15 days after delivery; for goods labeled with an expiration date, 15 days after use provided that such use occurs prior to the expiration date and for all other claims, 15 days after discovery of defect. Buyer shall afford Seller prompt and reasonable opportunity to inspect all goods as to which any claim is made. If Seller and Buyer are unable to reach settlement of any claim relating to goods covered hereby, Buyer shall institute legal action against Seller within one year after such claim arises and thereafter all such claims shall be barred notwithstanding any statutory period of limitations.

Indemnification

Buyer and Seller shall indemnify, hold harmless and defend each other from and against any and all actions, suits, liabilities, claims, fines, damages, losses and expenses including reasonable attorney's fees and expenses of litigation (the "Losses"), related to or arising from a) the injury of any person or the injury or destruction of any property arising out of and caused by the negligent acts or omissions of such party, any of its employees, subcontractors or authorized agents (collectively "Staff"); b) the material breach by such party of the terms and conditions of this Agreement; or c) violation of any law by such party or its Staff.

The indemnified party shall give the indemnifying party prompt notice of any claim that could give rise to a claim for indemnity under this Agreement and the indemnified party shall cooperate with the indemnifying party in the defense of any claim for which indemnity is provided. The indemnifying party shall be permitted to defend the claim and make all decisions thereto, including but not limited to hiring counsel of its choosing. The indemnifying party shall also have the sole right to settle any indemnified legal claim provided that it obtains a complete release for the indemnified party. This Section shall survive the termination or expiration of this Agreement for any reason.

Please see
<http://www.infraredx.com/infraredx-intellectual-property/>
for a listing of patents.

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Clinical Benefit

The clinical benefit provided by the Makoto™ Intravascular Imaging System is in the identification of the lipid core plaques (LCPs). By identifying LCPs, the Makoto™ Intravascular Imaging System provides information that is not as easily available from other technologies. The potential benefits of NIRS-IVUS-

guided management are: 1) determination of the length of artery to stent, 2) selection of intensity of medical management, 3) identification and mitigation of the risk of embolic peri-stenting infarction, and 4) identification of patients and plaques at increased risk of future major adverse cardiac events (MACE).

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Clinical Benefit

The clinical benefit provided by the Makoto™ Intravascular Imaging System is in the identification of the lipid core plaques (LCPs). By identifying LCPs, the Makoto™ Intravascular Imaging System provides information that is not as easily available from other technologies. The potential benefits of NIRS-IVUS-guided management are: 1) determination of the length of the artery to stent, 2) selection or intensity of medical management, 3) identification and mitigation of the risk of embolic peri-stenting infarction, and 4) identification of patients and plaques at increased risk of future major adverse cardiac events (MACE).

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Clinical Benefit

The clinical benefit provided by the Makoto™ Intravascular Imaging System is in the identification of the lipid core plaques (LCPs). By identifying LCPs, the Makoto™ Intravascular Imaging System provides information that is not as easily available from other technologies. The potential benefits of NIRS-IVUS-guided management are: 1) determination of the length of artery to stent, 2) selection of intensity of medical management, 3) identification and mitigation of the risk of embolic peri-stenting infarction, and 4) identification of patients and plaques at increased risk of future major adverse cardiac events (MACE).

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Clinical Benefit

The clinical benefit provided by the Makoto™ Intravascular Imaging System is in the identification of the lipid core plaques (LCPs). By identifying LCP, the Makoto™ Intravascular Imaging System provides information that is not as easily available from other technologies. The potential benefits of NIRS-IVUS-guided management are: 1) determination of the length of artery to stent, 2) selection or intensity of medical management, 3) identification and mitigation of the risk of embolic peri-stenting infarction, and 4) identification of patients and plaques at increased risk of future major adverse cardiac events (MACE).

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Clinical Benefit

The clinical benefit provided by the Makoto™ Intravascular Imaging System is in the identification of the lipid core plaques (LCPs). By identifying LCP, the Makoto™ Intravascular Imaging System provides information that is not as easily available from other technologies. The potential benefits of NIRS-IVUS-guided management are: 1) determination of the length of artery to stent, 2) selection of intensity of

medical management, 3) identification and mitigation of the risk of embolic peri-stenting infarction, and 4) identification of patients and plaques at increased risk of future major adverse cardiac events (MACE).

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NOTE:

Information related to the safety and clinical performance (SSCP) of the Makoto™ Intravascular Imaging System is available in the European database on medical devices (Eudamed), where it is linked to the Basic UDI-DI located beneath the bar code to (01) on the product label.

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1. Using sterile technique, open the pouch and remove catheter from packaging. DO NOT remove catheter from the hoop until the hydrophilic coating has been activated;

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Infraredx, Inc.