

NIPRO GM-30

Inflation Device

Key Benefits

Greater Versatility & Performance

- Maximum pressure reading: 30 ATM/441 PSI
- Flexible and kink resistant tubing
- Two tubing lengths available (31 cm and 76 cm) for a variety of interventional procedures
- 3-way and 4-way stopcock styles available
- Efficient pressurization to achieve desired inflation with fewer turns
- Dispersion volume of 25 mL allows for usage during various interventional procedures

Improved Ergonomics & Comfort

- Ergonomic handle designed for easy rotation and comfortable inflation and deflation
- Ergonomic lever design that fits comfortably in the hand to allow quick release with minimal effort
- Comfortable hand grip

Greater Visibility

- Fluorescent display allows for viewing in darkened conditions
- Easily readable numbers for volume and pressure readings



Get Comfortable. GM-30 Inflation Device.

Request Information at
info@nipromed.com | www.nipro.com | +1.908.393.7030



GM-30 Inflation Device

**Ergonomically-
designed
handle**



Less force and
fewer turns



Smooth inflation,
with 1:1 response*

*internal testing results on file

**Easy-to-use
release lever**




Minimal effort
required for deflation



Suitable for many balloon
sizes and diameters

**Options for
various
interventional
procedures**



Two stopcock styles:
3-way and 4-way



Two tubing lengths:
31 cm and 76 cm

The GM-30 inflation device features a 30 ATM pressure capacity with a 25 mL syringe volume (appropriate for usage with larger sized balloons).

Three-way and four-way stopcock positions and multiple tubing lengths (31 cm and 76 cm) allow for a variety of interventional procedures and techniques. **The GM-30 inflation device: The choice is yours. Get comfortable.**

Ordering Information

Product Name	Product Reference	Contents	GTIN
GM-30 Inflation Device	GM-30N	30atm-25cc Inflation Device. 31 cm, high pressure tubing and 3 position stopcock.	04543660022648
GM-30 Inflation Device	GM-30NL	30atm-25cc Inflation Device. 76 cm, high pressure tubing and 3 position stopcock.	04543660022662
GM-30 Inflation Device	GM-30NF	30atm-25cc Inflation Device. 31 cm, high pressure tubing and 4 position stopcock.	04543660022655
GM-30 Accessory Kit	GLGMN-YOK0A	A single box containing 1 sterile pouch with OKAY® II Y-Connector Set and 1 sterile pouch with a GM-30N Inflation Device.	04543660023508
GM-30 Accessory Kit	GLGMN-YOK0E	A single box containing 1 sterile pouch with OKAY® II Y-Connector Set w Ext. Tube and 1 sterile pouch with a GM-30N Inflation Device.	04543660023515
GM-30 Accessory Kit	GLGMF-YOK0A	A single box containing 1 sterile pouch with OKAY® II Y-Connector Set and 1 sterile pouch with a GM-30NF Inflation Device.	04543660023485
GM-30 Accessory Kit	GLGMF-YOK0E	A single box containing 1 sterile pouch with OKAY® II Y-Connector Set w Ext. Tube and 1 sterile pouch with a GM-30NF Inflation Device.	04543660023492

Distributed by



Nipro Medical Corporation
200 Crossing Boulevard
Bridgewater, NJ 08807
USA

T: +1.908.393.7030

Request Information at

info@nipromed.com
www.nipro.com

Imported by

Infraredx, a Nipro Company
28 Crosby Drive, Suite 100
Bedford, MA 01730
USA
T. +1.888.680.REDX (7339)
www.infraredx.com

©2022 Infraredx™, a Nipro Company
All rights reserved. MK0019rD