


Infraredx Clarispro™

is the only extended bandwidth catheter available bringing you a crisper image and greater depth of field than any other HD-IVUS catheter on the market. Sharper images allow for easier identification of lumen, determining presence of plaque, and overall better visual of the vessel structure.


The impressive pushability, deliverability, and extended working length of the Infraredx Clarispro™ HD-IVUS are ideal for peripheral use.




Key Features



Longest working length (160cm)
for reaching distant peripheral imaging sites



Extended pullback length (150mm)
for maximizing imaging in a single scan



Low crossing profile (3.2Fr)
for accessing hard-to-reach places like the tibial arteries




Fast pullback speed (up to 20x)
saving time and making each procedure more efficient




No additional accessories required
for greater savings and more efficient workflow



Best-in-class image clarity
to determine optimal plaque modification therapy




Clarispro IVUS Image
35-65 MHz extended bandwidth



Competitor IVUS Image
20 MHz

Infraredx Clarispro™ HD-IVUS

Technical Specifications



Physical Parameters

Working Length	Hydrophilic Coating	Tip-to-IVUS Image	Diameter of Catheter			Guiding Sheath/ Catheter Compatibility (Minimum Inside Diameter)	Optimal Guide Wire Size
			Entry Profile	Crossing Profile	Proximal Shaft		
160 cm	50 cm	16 mm	2.4 F (0.80 mm)	3.2 F (1.09 mm)	3.6 F (1.20 mm)	6F (0.067 in)	0.014 in (0.36 mm)

Imaging Parameters

IVUS Frequency Range	Maximum Image Diameter	Motorized Pullback	Motorized Live IVUS	Imaging Length	Indicated Vessels
35-65 MHz	16 mm	0.5, 1.0, 2.0 mm/s	0.0, 2.0, 10.0 mm/s bi-directional	150 mm	coronary and peripheral vessels, excluding cerebral vessels*

Ordering Information

Catalog Number	Description	GTIN
TVC-E195-42	Infraredx Clarispro™ 35-65 MHz HD-IVUS Catheter	00857595006180

**The Makoto Intravascular Imaging System™ is intended for ultrasound examination of coronary and peripheral intravascular pathology. Intravascular ultrasound imaging is indicated in patients who are candidates for transluminal coronary and peripheral interventional procedures. The System is not indicated for use in the cerebral vessels. See product instructions for full product description, contra-indications, warnings, and instructions.*

Increase the Capabilities of Your Makoto™ System

Infraredx Clarispro™ HD-IVUS works seamlessly with the Makoto Intravascular Imaging System – an imaging platform that benefits your whole practice to improve care for your patient’s polyvascular disease. The Makoto system is also compatible with the Dualpro IVUS+NIRS catheter to provide the only comprehensive coronary assessment tool FDA cleared to identify patients and lesions at increased risk of major adverse cardiac events.

Distributed by:

infraredx™
A NIPRO COMPANY

Infraredx, a Nipro Company
28 Crosby Drive, Suite 100
Bedford, MA 01730
T. +1.888.680.REDX (7339)

Request information at:
info@infraredx.com | www.infraredx.com

Imported by:
Infraredx, a Nipro Company
28 Crosby Drive, Suite 100
Bedford, MA 01730
USA
T. +1.888.680.REDX (7339)